

LDT 100

Operating Table



Product Description

- Adopt Taiwan imported TiMOTION electric push rod as the power source of the operation beds. The output power is large and quality is reliable
- Adopt German imported STABILUS as the power source of the operation beds. The output power is large and the quality is reliable.
- Operation beds adopt Medical Grade 304# full stainless steel materials, the surface is processed through matte processing, feeling more modern, more convenient to clean.
- Full set of micro electric control system with low noise, overload protection, safety reliability.
- The entire table of operating table can move forward and backward to realize full position C-arm and X-ray radiography.



LDT 100 Spesification

Head Board up	$\geq 45^{\circ}$ (up);	Dimension	2100x520mm
Head Board Down	$\geq 25^{\circ}$ (down)	Power Source	Electric
Back Board up	$\geq 75^{\circ}$ (up)	Voltage	110-130/220-240V
Back Board down	$\geq 15^{\circ}$ (down)	Material	3mm Stainless steel 304#
Leg Board Down	$\geq 90^{\circ}$	Lowest and Highest	680-980mm
Leg Board Outward	$\geq 90^{\circ}$	Elevation(up/down)	0-300mm
Kidney Bridge (optional)up/down	≥ 120 mm	Left/Right Lean	$\geq 20^{\circ}$
Max loading	320KGS	Forward Lean	$\geq 20^{\circ}$
		Backward Lean	$\geq 20^{\circ}$

Accessories	Quantity
Square Clamp I	5pcs
Square Clamp II	2pcs
Arm Rest	2pcs
Chest&Waist Supporter	2pcs
Leg-Crutches	2pcs
Shoulder Supporter	2pcs
Supporter	1pc
External Kidney Bridge	Optiona

